

New!

# ClearCalc

## Independent Calculation Software

Finalizing treatment plans is now easier than ever with ClearCalc, an independent secondary dose calculation software that delivers accurate results and clinical time savings. As a versatile solution, ClearCalc validates plan accuracy with direct plan import—eliminating the need to import/export plan data—for Eclipse users or by exporting DICOM files for those without Eclipse. Documentation is effortless as ClearCalc results can be included in ClearCheck plan reports with a single click.



- ✓ Full integration with ClearCheck streamlines planning workflows
- ✓ The automatic point selection workspace saves time and provides insightful pass/fail dose statistics
- ✓ Accepts DICOM imports from multiple TPS vendors for flexibility in mixed environments, including CyberKnife
- ✓ Supports 3D, IMRT, VMAT, SRS/SBRT, brachytherapy, electrons, & more in a single solution
- ✓ Diode support allows users to calculate expected doses for comparison against in-vivo measurements