

# ClearCalc Feature Compatibility

If you're unsure which ClearCalc features will be compatible with your clinic's current version of Eclipse, refer to this quick guide to find out.

ECLIPSE VERSION	11	13.5	13.6, 13.7	15.1, 15.5, 15.6
Auto commissioning	✓	✓	✓	✓
Auto-select point feature	✓	✓	✓	✓
DICOM compatible	✓	✓	✓	✓
Automatically send photon data to ClearCalc			✓	✓
Automatically send BrachyVision data to ClearCalc				✓
Automatically send electron data to ClearCalc				✓
Integration with ClearCheck (launch from ClearCheck and stream results back to ClearCheck report)			photon only	photon, brachytherapy, electrons