

# EZFluence Feature Compatibility

If you're unsure which EZFluence features will be compatible with your clinic's current version of Eclipse, refer to this quick guide to find out or ask one of our specialists.

Eclipse Version	11	13.5	13.6	13.7	15.1+
Automate field-in-field planning	✓	✓	✓	✓	✓
Automate e-comp planning	✓	✓	✓	✓	✓
DICOM compatible	✓	✓	✓	✓	✓
Automatically send plan data to EZFluence			✓	✓	✓
Support intermediate dose workflow between Eclipse and EZFluence doses			✓	✓	✓
Automatically send EZFluence data back to Eclipse					✓
Export auto-generated PTV_EVAL_EZ structure to Eclipse					✓
Automated IMRT-to-3D billing justification planning					✓